

Interim Report

Module-A. Disinfection efficacy of OptiMaser Microwave System.

Objective 4. Optimization of Corncob disinfection by OptiMaser.

Work Package- WP4.1-WP4.3

Introduction

Rodents should be housed with bedding because it allows foraging, burrowing, digging, and nest building and absorbs urine and faeces. Moreover, the type of bedding selected can influence animal wellbeing and experimental results. Bedding is evaluated for absorbency, biodegradability, toxicity, dust, palatability, comfort, cost, availability, damage to cage washers, and effect on research. Several bedding materials are available for rodents, including corncob, wood chips, paper products, and grass fibre pellets.

Untreated softwood bedding can affect rodent metabolism, and the aromatic hydrocarbons of cedar shavings can induce hepatic microsomal enzymes and cytotoxicity. Aspen bedding is associated with sneezing and respiratory pathology in rats. In a study comparing aspen shavings, virgin pulp, recycled paper, corncob, reclaimed wood pulp, virgin cellulose, pine shavings, and hardwood chips, ammonia was detectable last in corncob. Rats have been shown to prefer paper bedding to corncob; however, corncob has been used to minimize dampness and ammonia concentrations. Mice prefer material suitable for nest building such as cloth, cotton, or paper. Although corncob bedding is not a sterile material, it is a common bedding choice due to its absorbency, biodegradability, and ability to control ammonia levels

Corn cob advantages:

- Biodegradable
- Non-toxic
- Inert
- Absorbent-holding up to four times its weight of fluid
- Available in various uniform sizes
- Tough resilient & durable

Corn Cob bedding is also low dust, which is a big advantage to customers using individually Ventilated Cages (IVC's) and isolators, where HEPA filters often get clogged by conventional bedding. Corn Cob is also more densely packed, meaning that more cages can be filled from a single bag and the bags last

longer. Using Corn Cob bedding, rodent cages can be changed as infrequently as every 2 Week, as the cage remains dry and clean longer and labour is reduced as corn Cob bedding rarely gets wet or sticks to the cage requiring scrubbing or soaking. It is imperative to decontaminate animal house bedding before disposal. The animal house manure may contain high concern endocrine-disrupting compounds (EDCs) such as steroid hormones, rodents, such as mice and rats, can leave droppings that can spread bacteria, contaminate food, and cause allergic reactions. Microwave technology provides an efficient and quick solution to the animal bedding decontamination.

Material and Methods

A. Microorganisms used in the study

- Bacteria
 - *Staphylococcus aureus*
 - *Escherichia coli*
- Spores (*Bacillus subtilis*)
- Fungi
 - *Candida albicans*
 - *Candida glabrata*

B. Sample preparation

1. Used animal house Corncob bedding was procured from CSIR-IITR rodent animal house.
2. Corncob was packed into autoclavable biohazard plastic bags.
3. Water was added to these bags in 1:1 ratio (1kg Corncob/ 1lit water)
4. Bacterial culture were grown in LB broth.
5. Fungal culture were grown in YPD (Yeast Extract-1%, Dextrose-2% and Peptone- 2%).
6. *Bacillus subtilis* sporulation was performed in sporulation media (SM). Composition of SM used was- Nutrient broth- 0.8%, KCl- 0.1%, MgSO₄·7H₂O- 0.012%, pH adjusted to 7.6 with 1 M NaOH, volume adjusted to 1 litre with ddH₂O. Autoclaved and allowed to cool to 50°C. Just prior to use, following sterile solutions were added: Ca(NO₃)₂- 1mM, MnCl₂- 0.01mM, FeSO₄-0.001mM. 0.5 ml of an overnight grown *B subtilis* culture was inoculated in liquid DSM media in Erlenmeyer flasks. Flasks were kept at 37 °C in an orbital shaker at 200 rpm for 7 days. Generation of spores was monitored microscopically. After 7 days of incubation > 90% of spore population was observed.

7. All cultures were maintained at 37°C.
8. Viable cell density of at least 10^{10} cells was maintained for bacteria and fungus.

C. Methodology

1. 1g Corncob was taken in sterile 10 ml screw capped glass tubes and 1ml of microbial culture (10^{10} cells) was added to it. The microbial load was determined by total aerobic plate count or spectrophotometrically at 600 nm wavelength.
2. Tubes were deeply embedded into the Corncob bags and the bags were sealed with autoclave resistant rubber bands.
3. Bags were arranged in two rows within the OptiMaser bucket namely top and bottom.
4. Following OptiMaser treatment 100µl of this culture was plated in triplicates on suitable medium as untreated control.
5. For calculating the disinfection efficacy (Log reduction) total aerobic plate count procedure was followed.
6. For determining the Corncob disinfection efficacy, each organism was subjected to OptiMaser treatment within Corncob matrix at 100°C for two different time points (10 min and 30 min) (**Fig.1**).
7. For statistical significance disinfection efficacy of OptiMaser for each time point was determined using three biological replicates for each organism.
8. After OptiMaser treatment serial dilutions were prepared from each sample (up to 10^{-5}) and 100 µl sample of 10^{-1} , 10^{-3} , and 10^{-5} dilutions were plated in triplicates on suitable medium and incubated at 37°C for 48 h followed by calculation of log reduction (**Fig.1**).
9. As per 2016 guidelines of Biomedical waste management, *Bacillus atrophaeus* spores (10^4 spores) should be used as biological indicator for monitoring sterilization process. We used self-contained vials from gke for confirming the disinfection process in Corncob load.
10. Following OptiMaser cycles were used for this analysis:
 - Sprinkle time- 5 sec
 - Filament time- 90 sec
 - Microwave on time- 300- sec
 - Microwave Off time- 120 sec
 - Temp hold time- **600- 1800** sec
 - Cooling mode time- 400 sec

Water draining time- 0

Temp. Thresh hold- 100 °C

Printing numbers- 6

Results

Disinfection efficacy of Bacteria:

1. OptiMaser bacterial disinfection efficacy in Corncob matrix was tested using representative gram negative (*Escherichia coli* K12) and gram positive bacteria (*Staphylococcus aureus*).
2. Bacteria was exposed to OptiMaser at 100°C for 10 and 30 min respectively.
3. It was observed that at 10 min OptiMaser treatment time the disinfection of bacteria was good with 10 log reduction achieved in the top row. The log reduction was very reasonable in the bottom row where 6 log reduction was achieved for both gram negative and gram positive bacteria (Fig. 2).
4. After 30min or 40 min of OptiMaser exposure time a 10 log reduction was achieved in both top and bottom row for both gram negative and gram positive bacteria (Fig. 3 & 4).
5. Total aerobic plate count method was used for determine disinfection efficacy

Disinfection efficacy of Fungi:

1. OptiMaser Fungal disinfection efficacy in Corncob matrix was tested using representative pathogenic fungi *Candida albicans* and *Candida glabrata*.
2. It was observed that at 10 min OptiMaser treatment time the disinfection of Fungi was good with 10 log reduction in top row. However the log reduction was reasonable in the bottom row where 6 log reduction was achieved for both *C. albicans* and *C. glabrata* (Fig. 2).
3. After 30 min or 40min of OptiMaser exposure time a 10 log reduction was achieved in both top and bottom row for both *C. albicans* and *C. glabrata* (Fig. 3 & 4).
4. Total aerobic plate count method was used for determine disinfection efficacy

Disinfection efficacy of Spores:

1. OptiMaser heat resistant bacteria spore disinfection efficacy in Corncob matrix was tested using representative *Bacillus subtilis* spores.

2. It was observed that at 10 min OptiMaser treatment time the disinfection of bacterial spores was very poor with only 2 log reduction in top row. However the spores were not killed in the bottom row (Fig. 2).
3. After 30 min of OptiMaser exposure time a 10 log reduction of bacterial spores was achieved in the top row whereas 4 log reduction was observed in the bottom row (Fig. 3).
4. After 40 min of OptiMaser exposure time a 10 log reduction of bacterial spores was achieved in the top row whereas 6 log reduction was observed in the bottom row (Fig. 4).
5. Total aerobic plate count method was used for determine disinfection efficacy.
6. OptiMaser exposure at 100°C for 40 min was good enough to kill 10^6 *B. atrophaeus* spores, as observed by no change in colour of the growth media after 48hrs of incubation at 37°C. The colour of growth medium of untreated vial turned yellow within 8hrs of incubation at 37°C (Fig. 4D).

Conclusions

At least 10 log disinfection efficacy of representative bacterial and fungal culture in Corncob matrix was achieved via OptiMaser treatment at 100°C with a hold time of 40 min. A disinfection efficacy of 10 log reduction was achieved for heat resistant bacterial spores in the top row. Acceptable disinfection efficacy of bacterial spores (6 log reduction) was observed in the bottom row. Self-contained biological indicator having 10^6 *B. atrophaeus* spores was used for monitoring disinfection process. The spore vial was subjected to 40min OptiMaser exposure at 100°C along with Corncob load. No growth was recorded in the spore vial after 48hrs incubation at 37°C, suggesting successful disinfection process.





Figure 1: Plan of study for determining the Corncob disinfection efficacy of OptiMaser-30 (2.45 GHz and 1.5 kW). (A) Disinfection efficacy of OptiMaser was determined in Corncob matrix using representative bacteria, fungi and spores as indicated in the figure. (B) Due to complex matrix microbe disinfection efficacy was calculated at 100°C (two hold times). (C) For statistical significance three biological replicates of each samples were exposed to OptiMaser, following microwave exposure, samples were serially diluted and 10^{-1} , 10^{-3} and 10^{-5} dilutions were plated in triplicates on LB agar plates for total aerobic plate. (D) Brief methodology. Bacteria was grown in LB media at 37°C and 10^{10} cells were exposed to OptiMaser. Fungus was grown in YPD (Yeast Extract-1%, Dextrose-2%, Peptone- 2%) broth at 37°C and 10^{10} cells were exposed to OptiMaser. *Bacillus subtilis* spores were generated via magnesium sulphate method.



B

| Temp. hold mode | Max . Temp. achieved | Total run time | Log reduction | | | | | |
|------------------|----------------------|----------------|---------------|--------|-------|--------|--------|--------|
| | | | Bacteria | | Fungi | | Spores | |
| | | | Top | Bottom | Top | Bottom | Top | Bottom |
| 100°C for 10 min | 100°C | 32min | 10 | 6 | 10 | 6 | 2 | - |

C

| Log Reduction | % Reduction of bacteria |
|---------------|-------------------------|
| 10 | 99.99999999 |





Figure 2. Corncob Disinfection efficacy of OptiMaser at 100°C. (A) Representative plates for total aerobic plate count. **(B)** OptiMaser parameter at 10min hold time at 100°C and Log reduction in bacteria, spores and Fungi upon OptiMaser treatment. **(C)** Relation between log reduction and % reduction. **(D)** Arrangement of husk bags in OptiMaser bucket. Bags were arranged in two rows namely top and bottom.



B

| Temp. hold mode | Max . Temp. achieved | Total run time | Log reduction | | | | | |
|------------------|----------------------|----------------|---------------|--------|-------|--------|--------|--------|
| | | | Bacteria | | Fungi | | Spores | |
| | | | Top | Bottom | Top | Bottom | Top | Bottom |
| 100°C for 30 min | 100°C | 53min | 10 | 9 | 10 | 10 | 10 | 4 |

C

| Log Reduction | % Reduction of bacteria |
|---------------|-------------------------|
| 10 | 99.99999999 |





Figure 3. Corncob Disinfection efficacy of OptiMaser at 100°C. (A) Representative plates for total aerobic plate count. (B) OptiMaser parameter at 30min hold time at 100°C and Log reduction in bacteria, spores and Fungi upon OptiMaser treatment. (C) Relation between log reduction and % reduction. (D) Arrangement of husk bags in OptiMaser bucket. Bags were arranged in two rows namely top and bottom.



B

| Temp. hold mode | Max . Temp. achieved | Total run time | Log reduction | | | | | |
|------------------|----------------------|----------------|---------------|--------|-------|--------|--------|--------|
| | | | Bacteria | | Fungi | | Spores | |
| | | | Top | Bottom | Top | Bottom | Top | Bottom |
| 100°C for 40 min | 100°C | 65 | 10 | 10 | 10 | 10 | 10 | 6 |

C

| Log Reduction | % Reduction of bacteria |
|---------------|-------------------------|
| 10 | 99.99999999 |



Figure 4. Corncob Disinfection efficacy of OptiMaser at 100°C. (A) Representative plates for total aerobic plate count. (B) OptiMaser parameter at 40min hold time at 100°C and Log reduction in bacteria, spores and Fungi upon OptiMaser treatment. (C) Relation between log reduction and % reduction. (D) Disinfection of 10^6 *B. atrophaeus* spores upon OptiMaser exposure at 100°C for 30min along with rice husk load. No change in colour of the medium after 48hrs of incubation at 37°C suggest complete disinfection. C, Untreated control; T, Treated vial.