



Specifically formulated for

## HOSPITAL & HEALTH CARE

Disinfectant | Sanitizer | Bactericide | Broad Spectrum Virucide

**CRITICAL DISINFECTION:** Use EnviroNize HAISAN™ to disinfect areas such as medical equipment, accessories, patient rooms, reception areas, high-touch surfaces, bathroom fixtures, and any washable food and non-food contact surface where disinfection is required.

### PURE HYPOCHLOROUS

pH | PPM OPTIMIZED  
DIN 02527553



**KILLS 99.999%**

✓ Bacteria ✓ Viruses ✓ Germs

DEODORIZER

FRAGRANCE FREE

#### EFFECTIVE AGAINST

- Human Coronavirus
- Human Adenovirus type 5
- Staphylococcus aureus
- Pseudomonas aeruginosa
- Enveloped and Non-Enveloped viruses
- Other Bacteria

# Proactive Prevention of Pathogens™

## Sustainable Disinfectant

EnviroNize HAISAN™ brand of HOCL (Hypochlorous Acid) is uniquely certified organic input, single ingredient, effective, an environmentally alternative disinfectant and Made in Canada.


- Non-WHMIS regulated
- Non-Hazardous
- Zero Surface Contact Residue



Health	0 - None
Flammability	0 - None
Physical	0 - None
Personal Protection	A - Splash Goggles

## Powerful Disinfectant

EnviroNize® HAISAN™ (Hypochlorous Acid) can penetrate pathogenic microorganism cell walls very easily, making it a very effective disinfectant.



## RTU (Ready to Use) Disinfectant

EnviroNize® HAISAN™ is a Hospital/Healthcare Disinfectant/Sanitizer, Bactericide and a Broad Spectrum Viricide

EnviroNize® HAISAN™ is effective against, but not limited to, Human Coronavirus, Human Adenovirus type 5, Staphylococcus aureus, Pseudomonas aeruginosa and enveloped and non-enveloped viruses.

EnviroNize® HAISAN™ is to be used as the primary disinfectant for hard, non-porous, non-food contact surfaces in healthcare practices, including ACUTE CARE, ALTERNATE CARE and LONGTERM CARE areas.

EnviroNize® HAISAN™ may be used to pre-clean or decontaminate critical or semi critical devices prior to sterilization or high-level disinfection.

EnviroNize® HAISAN™ is for use on hard non-porous, non-food contact inanimate surfaces. EnviroNize® HAISAN™'s predominant antimicrobial agent is hypochlorous acid, an efficient and efficacious species of chlorine. Hypochlorous acid kills bacteria and viruses. EnviroNize® HAISAN™ is an oxidized pH neutral water-based solution that is Ready to Use (RTU).

# What Makes it the Best Disinfectant Choice?

EnviroNize™ **HAISAN™** DIN 02527553  
Re-Imagining Clean

- IS A CONCENTRATED ELECTRO-CHEMICALLY ACTIVATED SOLUTION rich in oxygen, in which pathogens cannot survive.
- IS A POWERFUL ANTIMICROBIAL SOLUTION that is both biodegradable and non-accumulative. The key active ingredient, hypochlorous acid, is a naturally occurring molecule synthesized from an electrolyzed solution of salt and water.
- IS NON-WHIMS REGULATED at pH Neutral
- CAN PENETRATE PATHOGENIC MICROORGANISM CELL WALLS
- IS HIGH IN ORP SOLUTION (ORP above +900mV)
- ATTACKS AN ARRAY OF MICROBIAL ORGANISMS, a process that is chemically selective with no toxic or mutagenic side effects for higher forms of life.
- IS SAFE WHEN USED AROUND PLANTS, ANIMALS, MARINE OR HUMANS
- IS AN ENVIRONMENTALLY CONSCIOUS AND SUSTAINABLE SOLUTION



**ALL HARD NON-POROUS SURFACES**



**PROACTIVE  
PREVENTION  
OF PATHOGENS™**

**ENVIRONMENTAL**

EnviroNize HAISAN™ is made from all natural ingredients, biodegradable and 100% eco friendly.

**SOCIAL**

EnviroNize HAISAN™ is non-irritant to the skin, non-flammable, does not promote skin dryness, and is viewed as a positive ecological impact.

**ECONOMIC**

EnviroNize HAISAN™ provides refillable and re-useable containers and devices significantly reducing your carbon footprint.

## Powerful & Sustainable HOCl Technology



**HOSPITAL-  
ACQUIRED  
INFECTION**  
& its prevention

Electrolyzed sodium hypochlorous acid (HOCl), produces an efficient disinfecting adjunct to standard terminal cleaning and provides many advantages in optimizing infection prevention in the hospital environment.

**SOURCE:** [www.infectioncontroltoday.com/view/hypochlorous-acid-definitive-terminal-cleaning-hospital-environment](http://www.infectioncontroltoday.com/view/hypochlorous-acid-definitive-terminal-cleaning-hospital-environment)



# EnviroNize® HAISAN™ Accessories

## DEPLOY WITH THE FOLLOWING DEVICES:

- EnviroNizer Airless Pump Sprayer
- EnviroNize ULV
- EnviroNizer Cordless
- EnviroNize Ultra Max
- EnviroNize Vapore

- Trigger-based dispensing technology
- Continuous spray
- Re-useable head and convenient refill with EnviroNize® HAISAN™
- Disperse | Refill | Re-use



**ITEM # HAIS5007**  
EnviroNize® HAISAN™  
185 ml Airless Pump Sprayer

## ECONOMICAL

Refill & re-use from a larger EnviroNize® HAISAN™ container



## Micro Fibre Towels

- Prevents cross contamination
- Separates cleaning & disinfection



**ITEM # ENVFT-A**   **ITEM # ENVFT-C**  
Anolyte                      Catholyte



**ITEM # ENVSPIGOT**  
3/4" Natural Pour Spigot

## ITEM # ENV3SHELF

Store & dispense bulk quantities of EnviroNize® solutions



# EnviroNize® HAISAN™ Solutions

ITEM	DESCRIPTION	CASE WEIGHT	UPC CODE	PACK SIZE	CASES FOR PALLET
EHAS5007	185 ml   6.25 oz Haisan Disinfectant w/EnviroNizer Airless Pump Sprayer	4.15 lbs   1.88 kg	625570900736	6	77 (7-11)
EHAS5002	1000 ml   35.20 oz Haisan Disinfectant	16.0 lbs   7.26 kg	625570100891	6	154 (7-22)
EHAS55002-TS	1000 ml   35.20 oz Haisan Disinfectant w/Trigger Sprayer	16.0 lbs   7.26 kg	625570900682	6	154 (7-22)
EHAS5003	3785 ml   128 oz Haisan Disinfectant	36.2 lbs   16.42 kg	625570100884	4	72 (6-12)
EHAS5004	20 L Haisan Disinfectant	47.00 lbs   19.62 kg	625570100877	1	48 (4-12)



Made in Canada, Manufactured under License from PJ Holdings Corp.  
11-1673 Richmond St. PMB 110, London, ON, Canada N6G 2N3  
877.432.8548 | [haisan@environize.ca](mailto:haisan@environize.ca)  
[www.environize.ca/haisan](http://www.environize.ca/haisan)

